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Shunt Band

### 1. Device Identification

Device Name: Shunt Band  
 Medical Device Class: Class I  
 Intended Users: Healthcare professionals only  
 Reusable: Yes

## 2. Intended Purpose

The device is intended for the temporary fixation of a swab (gauze pad) on a puncture site to apply defined pressure to the skin, thereby supporting hemostasis following a puncture.

It is used as an aid in procedures where manual compression is required for bleeding control, particularly in the context of dialysis treatment or other venous punctures.

The device is intended for use on intact skin and does not penetrate the body.

## 3. Intended User and Patient Population

### Intended User:

The device is intended solely for use by trained healthcare professionals who are familiar with the proper application of compression aids and/or shunt management.

### Intended Patient Population:

The device may be used on individuals of all age groups for whom manual compression of a puncture site is clinically indicated.

## 4. Duration and Frequency of Use

Contact Duration: Short-term use for hemostasis, typically 2-5 minutes, at most until secure hemostasis is achieved.

Frequency: Based on clinical judgment and the number of punctures performed.

Body Contact: The device must only come into contact with intact skin; the puncture site itself must be completely covered with a swab.

## 5. Contraindications

The device must not be used:

- On inflamed, infected, or severely damaged skin.
- On patients with known skin hypersensitivity to the materials of the device.
- In the presence of signs of relevant circulatory disorders.

## 6. Warnings and Precautions

- The device is intended for professional use only.
- Ensure that the puncture site is completely covered with a swab before the device is applied.
- Important: While the device must generate sufficient pressure for hemostasis, it must not completely occlude the shunt or the blood flow. Check the shunt signal after application.
- Inspect the device before each use. Do not use if it is damaged, soiled, or no longer functional.
- Do not use the device for purposes other than fixation of a swab on a puncture site as described in this IFU.
- Do not modify the device.
- Users must ensure that the current version of the electronic instructions for use is consulted prior to use of the device.

## 7. Residual Risks

Residual risks associated with the use of the device are considered low and mainly relate to:

- Inappropriate compression (too tight, which can promote shunt thrombosis, or too loose, leading to hematomas).
- Skin irritation in the contact area of the device.

These risks are mitigated by:

- Use by trained healthcare professionals only
- Adherence to the instructions and warnings provided in this IFU

When used as intended, the benefits of the device outweigh the remaining residual risks.

## 8. Instructions for Use

- Place a swab completely over the puncture site.
- Apply the device such that it securely fixes the swab and exerts even, light pressure.
- Check after 2-5 minutes whether hemostasis has been achieved.
- Remove the device as soon as there is no further bleeding.

## 9. Cleaning and Maintenance

### Cleaning

Clean the device after each use with a mild detergent and disinfectant and water at max. 60 °C (manual wash) or a disinfectant suitable for this medical device.

Ensure the device is completely dry before reuse.

## Maintenance

- No maintenance is required.

## 10. Storage and Handling

- Store at room temperature in a clean, dry environment.
- Protect from excessive heat, moisture, and mechanical stress.
- Do not store in direct sunlight.

## 11. Electronic Instructions for Use (eIFU)

The instructions for use for this device are provided exclusively in electronic form.

The current electronic IFU is available at:

<http://doc.1m4.eu/a4.pdf>

A paper version of the IFU is available upon request for healthcare professionals.

## 12. Disposal

Dispose of the device in accordance with local regulations for medical devices.

## 13. Traceability

Each device is identified by a batch (LOT) number. The batch number must be documented to ensure traceability in accordance with applicable regulatory requirements.

## 14. Reporting of Incidents

Any serious incident related to the use of this device must be reported to the manufacturer and the competent authority in accordance with applicable regulations.

## 15. Explanation of Symbols

The symbols on the label are listed in the table at the end of this IFU.

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